
| | | |
|-----------------------------------|--|----|
| Instructions for use | HAMILTON-BC8010, Breathing circuit set | 6 |
| Gebrauchsanweisung | HAMILTON-BC8010, Beatmungsschlauchset | 10 |
| Instructions d'utilisation | Kit de circuit respiratoire HAMILTON-BC8010 | 14 |
| Instrucciones de uso | del equipo respiratorio HAMILTON-BC8010 | 18 |
| Istruzioni per l'uso | HAMILTON-BC8010, set circuito paziente | 22 |
| Instruções de uso | para HAMILTON-BC8010, kit de respiração | 26 |
| 使用说明 | HAMILTON-BC8010, 呼吸装置 | 30 |
| 使用説明書 | HAMILTON-BC8010、呼吸回路セット | 33 |
| Инструкция по эксплуатации | комплекта дыхательного контура HAMILTON-BC8010 | 37 |
| Kullanım Talimatları | HAMILTON-BC8010, Solunum devresi seti | 42 |

REF 260185


CE 0197

ERC

HAMILTON
MEDICAL



Breathing circuit diagram | Schematische Darstellung des Beatmungsschlauchsystems | Schéma du circuit respiratoire | Diagrama del circuito respiratorio | Diagramma del circuito paziente | Esquema do circuito de respiração | 呼吸管路图表 | 呼吸回路の図 | Схема дыхательного контура | Solunum devresi şeması

- 1** Ventilator inspiratory limb
Inspirationsschenkel für das Beatmungsgerät
Branche inspiratoire du ventilateur
Rama inspiratoria del respirador
Branca inspiratoria del ventilatore
Alça inspiratória do respirador
呼吸机吸气肢
人工呼吸器側の吸気回路
Патрубок вдоха со стороны аппарата ИВЛ
Ventilator inspiratuvur parçası
- 2** Humidifier chamber
Befeuchterkammer
Réservoir de l'humidificateur
Cámara del humidificador
Camera dell'umidificatore
Câmara do umidificador
湿化器水箱
加温加湿器チャンバ
Камера увлажнителя
Nemlendirici haznesi
- 3** Patient inspiratory limb, heated
Inspirationsschenkel zum Patienten, beheizt
Branche inspiratoire côté patient, chauffante
Rama inspiratoria del paciente, con calefacción
Branca inspiratoria del circuito paciente, riscaldada
Alça inspiratória do paciente, aquecida
病人吸气肢, 加热
患者側の吸気回路、加熱
Патрубок вдоха со стороны пациента (нагреваемый)
Hasta inspiratuvur parça, ısıtılmı
- 4** Connector: OD10/OD10
Anschluss: AD10/AD10
Connecteur: DE 10/DE 10
Conector: DE10/DE10
Connettore: D.E. 10/D.E. 10
Conector: DE10/DE10
接头: 外径 10/ 外径 10
コネクタ: 外径 10/ 外径 10
Коннектор: НД: 10 / НД: 10
Konektör: DÇ 10/DÇ 10
- 5** Unheated inspiratory limb extension, for use in incubator
Nicht beheizte Verlängerung für den Inspirationsschenkel, zur Verwendung im Inkubator
Extension de branche inspiratoire non chauffante, à utiliser en couveuse
Ampliación de rama inspiratoria sin calefacción para uso en incubadora
Prolunga della branca inspiratoria non riscaldata, per l'uso dentro l'incubatrice
Extensão de alça inspiratória não aquecida para ser usada na incubadora
在保温箱内使用未加热的延长吸气肢
加熱されていない吸気エクステンションチューブ、インキュベーター内で使用
Удлинитель ненагреваемого патрубкa вдоха, для использования в кувезах
Isitmali olmayan inspiratuvur parça uzantısı; inkübatörde kullanılmı içindir.
- 6** Y-piece adapter with sampling port OD10/OD10/D15
Y-Stück-Adapter mit Probenentnahmeanschluss AD10/AD10/D15
Adaptateur pour pièce en Y avec port d'échantillonnage DE 10/DE 10/DI 15
Adaptador de pieza en Y con puerto de muestreo DE10/DE10/DI15
Adattatore del raccordo a Y con porta di campionamento D.E. 10/D.E. 10/D.I. 15
Adaptador peça em "Y" com porta de amostras DE10/DE10/DI15
Y 形管接头, 带采样端口外径 10/ 外径 10/ 内径 15
サンプリングポート付き Y ピースアダプタ、外径 10/ 外径 10/ 内径 15
У-образный коннектор с пробоотборным отверстием: НД: 10 / НД: 10 / ВД: 15
Y parçası adaptörü; DÇ 10/DÇ 10/Ç 15 örnekleme portu ile birlikte

- 7** Patient expiratory limb
 Expirationsschenkel für den Patienten
 Branche expiratoire patient
 Rama espiratoria del paciente
 Branca espiratoria del circuito paziente
 Alça expiratória do paciente
 病人呼气肢
 患者側の呼吸回路
Патрубок выдоха со стороны пациента
 Hasta ekspiratuvar parçası
- 10** Luer connector, male, symmetrical (2x)
 Luer-Stecker, symmetrisch (2x)
 Connecteur Luer, mâle, symétrique (2x)
 Conector Luer, macho, simétrico (2x)
 Connettore Luer, maschio, simmetrico (2x)
 Conector luer, macho, simétrico (2x)
鲁尔接头，外螺纹，对称型 (2x)
 ルアーコネクタ、オス、対称型 (2個)
Штекерный люэровский наконечник, симметричный (2x)
 Luer konektör, erkek, simetrik (2x)
- 8** T-piece OD10/ID10 with pressure line, 1.6 m
 T-Stück AD10/ID10 mit Druckmessleitung, 1,6 m
 Pièce en T DE 10/DI 10 avec ligne de pression de 1,6 m
 Pieza en T DE10/DI10 con línea de presión, 1,6 m
 Raccordo a T, D.E. 10/D.I. 10, con linea di pressione, 1,6 m
 Peça em "T" DE10/DI10 com linha de pressão, 1,6 m
T形管外径 10/ 内径 10, 带有耐压管路 (长 1.6 米)
 Tピース、外径 10/ 内径 10、圧力ライン付き、1.6 m
T-образный коннектор, НД: 10 / ВД: 10, с линией контроля давления 1,6 м
 T parçası DÇ10/İÇ10; 1,6 m'lik basınç hattı ile birlikte
- 11** Soft adapter: ID15/ID10 (2x)
 Weicher Adapter: ID15/ID10 (2x)
 Adaptateur souple : DI 15/DI 10 (2x)
 Adaptador blando : DI15/DI10 (2x)
 Adattatore morbido: D.I. 15/D.I. 10 (2x)
 Adaptador macio: DI15/DI10 (2x)
柔性接头：内径 15/ 内径 10 (2x)
 ソフトアダプタ：内径 15/ 内径 10 (2 個)
Мягкий адаптер: ВД: 15 / ВД: 10 (2x)
 Yumuşak adaptör: İÇ15/İÇ10 (2x)
- 9** Elbow adapter, 90°: ID22/ID15/OD22
 Ellenbogenadapter, 90°: ID22/ID15/AD22
 Adaptateur coudé, 90°: DI 22/DI 15/DE 22
 Adaptador acodado, 90°: DI22/DI15/DE22
 Adattatore a gomito, 90°: D.I. 22/D.I. 15/D.E. 22
 Adaptador de cotovelo, 90°: DI22/DI15/DE22
弯接头，90°：内径 22/ 内径 15/ 外径 22
 エルボアダプタ、90°：内径 22/ 内径 15/ 外径 22
Угловой адаптер, 90°: ВД: 22 / ВД: 15 / НД: 22
 Dirsek adaptörü, 90°: İÇ22/İÇ15/DÇ22
- 12** Straight connector: OD10/ID15 (2x)
 Gerader Anschluss: AD10/ID15 (2x)
 Connecteur droit : DE 10/DI 15 (2x)
 Conector plano: DE10/DI15 (2x)
 Connettore diritto: D.E. 10/D.I. 15 (2x)
 Conector reto: DE10/DI15 (2x)
直套管：外径 10/ 内径 15 (2x)
 ストレートコネクタ：外径 10/ 内径 15 (2 個)
Прямой коннектор: НД: 10 / ВД: 15 (2x)
 Düz konektör: DÇ10/İÇ15 (2x)

Incubator and warming bed use | Verwendung eines Inkubators und eines Wärmebetts | Utilisation de la couveuse et du lit chauffant | Uso del incubador y la cama calefactora | Uso nell'incubatrice e nel letto riscaldante | Uso da incubadora e cama de aquecimento | 保温箱和保暖床的使用 | インキュベーターおよび温熱ベッドでの使用 | Использование кувеза и согревающей кровати | İnkübatör ve ısıtıcı yatak kullanımı

Conventions

⚠ WARNING

Alerts the user to the possibility of injury, death, or other serious adverse reactions associated with the use or misuse of the device.

⚠ CAUTION

Alerts the user to the possibility of a problem with the device associated with its use or misuse, such as device malfunction, device failure, damage to the device, or damage to other property.

NOTICE

Emphasizes information of particular importance.

1 Overview

The HAMILTON-BC8010 breathing circuit set (PN 260185) can only be used with the HAMILTON-H900 humidifier. For an overview of the breathing circuit set and its components, see page 2.

1.1 Intended use

The HAMILTON-BC8010 breathing set is intended to be used together with compatible Hamilton Medical respiratory gas humidifiers during invasive and noninvasive mechanical ventilation of pediatric and neonatal patients.

2 Safety information

⚠ WARNING

- Reprocessing of Hamilton Medical single-use products can affect the product properties and may cause injury to the patient. For example, a change to the surface structure during reprocessing may lead to a change in the tear strength or cause actual cracking. Furthermore, an altered surface structure

may result in, for example, a microbial aggregation of spores, allergens, and pyrogens, or cause an increase in the number of particles released as a result of chemical changes in the material properties.

- Hamilton Medical does *not* assume any warranty for the proper function of single-use items if they are reprocessed and reused by the user. Any liability on the part of Hamilton Medical for any damages resulting from the use of reprocessed or reused single-use parts, or use of any non-Hamilton Medical accessories and consumables, is excluded.
- The breathing circuit set must *not* be covered by any objects, such as sheets, towels, and so on.
- Use caution when positioning the breathing circuit. Place the breathing circuit in such a way that condensation cannot reach the patient.
- Follow hospital infection control procedures for all tests and use.

- Discard the product if it fails two subsequent preoperational checks or if there is any sign of damage.
 - Adding attachments or other parts/assemblies to a breathing system can change the pressure gradient across the ventilator, which can adversely affect ventilator performance and patient safety.
 - Connect only a single breathing circuit at a time. Do *not* connect any circuit as an extension to another one. Doing so can seriously harm or kill the patient.
 - Route the breathing circuit without tension and without any kinks from the ventilator or humidifier to the patient, and protect against unintended influences.
 - Do *not* use heated breathing circuits without existing gas flow. If the gas supply is interrupted, the humidifier must be turned off or be in **Standby** mode.
 - Do *not* touch the hot plate or the bottom of the chamber. The surfaces can reach a temperature of over 85°C. These hot surfaces radiate heat.
 - Heated limbs must *not* be placed directly on the patient's skin. Attach breathing circuits or tubing holders appropriately to avoid mechanical forces on the ET tube.
 - Only fill the humidifier chamber with sterile, demineralized water that meets the hospital's hygiene requirements.
 - Do *not* operate the device if the water level in the humidifier chamber exceeds the marked maximum.
 - Check the breathing circuit for condensation regularly and drain if required.
 - The humidifier must always be positioned below patient level. Do *not* operate the humidifier at an angle in excess of 10°.
 - Do *not* use this breathing circuit set with humidifiers other than HAMILTON-H900.
 - During high-frequency oscillation (HFO), regularly check the patient condition and ventilation settings, and refer to the HFO device Operator's Manual.
 - *The temperature symbol on the temperature probe must always face up and be visible.*
 - *Use the unheated limb extension only inside the incubator.*
 - *Always place the heated limb with the temperature probe outside the incubator.*
-
- NOTICE**
- Before use, read the *HAMILTON-H900 Humidifier Instructions for Use*.
 - The refill water should *not* be warmer than 37°C.
 - Follow the *Instructions for use* of the corresponding ventilator when using a pressure line.
 - Ensure that the water supply to the humidifier chamber is functioning properly.
 - Set the appropriate ventilator alarms.
 - Do *not* apply pressure to the water reservoir once the water refill has automatically stopped and the water has reached the maximum level.
 - Ensure the humidifier is placed in **Standby** or turned off before disconnecting any components.
 - Do *not* add any drugs or medication directly to the water in the water chamber. If the HAMILTON-H900 is used in combination with any medical

⚠ CAUTION

- *(USA only) Federal law restricts this device to sale by or on the order of a physician.*
- *Check the breathing set for damage prior to use. Discard the breathing set if there is any sign of damage.*
- *Handle used products as contaminated goods according to local laws and regulations or hospital internal procedures.*

gases or nebulized medications, follow the Instructions for use of the supplied application and make sure that it is suitable for use with active humidification.

- Use the stickers provided to indicate the date on which to replace the breathing circuit set.
- Check the stability of all connections, prior to use.
- Incompatible parts can result in degraded performance, which can affect safety.

3 About the breathing circuit and breathing circuit sets

The breathing circuit comprises the limbs with integrated heater wire and integrated temperature probe. The breathing circuit set comprises the breathing circuit and the humidifier chamber with water feed line.

4 Connecting the breathing circuit

Note that the limb connectors on the humidifier combine electrical connections with breathing circuit connectors. Ensure proper orientation of electrical contacts on the breathing circuit connectors to match the connecting element on the humidifier. For an overview of the breathing circuit set and its components, see page 2.

To connect the breathing circuit set to the humidifier

- 1 Insert the humidifier chamber completely into the humidifier until it clicks into place.
To remove it, pull the humidifier chamber out of the humidifier.
- 2 Connect the ventilator inspiratory limb first to the humidifier, and then to the inspiratory port on the ventilator.
- 3 Connect the ventilator expiratory limb to the expiratory port on the ventilator and the Y-piece.
- 4 Connect the patient inspiratory limb first to the humidifier chamber, and then to the Y-piece.
- 5 Insert the humidifier chamber water feed spike into the water reservoir.

5 Delivery

The components are delivered clean and ready for clinical use.

6 Disposal

A used component must be handled as contaminated. Follow all local, state, and federal regulations with respect to waste management and environmental protection when disposing of used parts.

7 Technical specifications

| Parameter | Specification ¹ |
|--|----------------------------|
| Flow range | 1 to 30 l/min |
| Connectors ^{2,3} Ventilator: Patient: | OD15 ID15 |
| Adapters Elbow adapter, 90°: | ID22/ID15/OD22 |
| Y-piece adapter with sampling port: | OD10/OD10/ID15 |
| Straight connector (2x): | OD15/OD10 |
| T-piece with pressure line 1.6 m: | OD10/ID10 |
| Luer connector, male, symmetrical (2x) | |
| Compliance (at 60 cmH ₂ O) | 1.3 ml/cmH ₂ O |
| Compressible volume | 800 ml |
| Smallest internal diameter | 8.00 mm |
| Water volume (min/max) | 80 ml / 200 ml |

| Parameter | Specification ¹ |
|---|--|
| Flow resistance (at 2.5 l/min) | |
| Inspiratory: | 0.03 cmH ₂ O/(l/min) |
| Expiratory: | 0.02 cmH ₂ O/(l/min) |
| Flow resistance (at 15 l/min) | |
| Inspiratory: | 0.06 cmH ₂ O/(l/min) |
| Expiratory: | 0.04 cmH ₂ O/(l/min) |
| Gas leakage (at 60 cmH ₂ O) | < 30 ml/min |
| Maximum operating pressure | 20 kPa |
| Circuit length | 1.65 m (approx. 5.5 ft) |
| Compatibility | HAMILTON-H900 humidifier |
| Operating tempera- ture | Recommended: 18°C to 26°C (65°F to 78°F) |
| Storage temperature | -20°C to 60°C (-4°F to 140°F) |
| Storage and operat- ing humidity | 5% to 95% relative humidity, noncon- densing |
| Maximum operating duration: | 28 days |

| Parameter | Specification ¹ |
|-------------|---|
| Conforms to | IEC 60601-1:2005/ A1:2012 ISO 80601-2- 74:2017 EN ISO 5356-1:2015 EN ISO 5367:2014 IEC 62366-1:2015 ISO 10993-1:2018 ISO 13485:2016 |

1. All values presented are worst case scenarios
2. Internal diameter (ID) and outer diameter (OD) in mm.
3. 15 mm and 22 mm conical connectors according to ISO 5356-1.

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For devices manufactured in Switzerland

Manufacturer

Hamilton Medical AG
Via Crusch 8, 7402 Bonaduz, Switzerland
☎ +41 58 610 10 20
info@hamilton-medical.com
www.hamilton-medical.com

medin Medical Innovations GmbH
Adam-Geisler-Str. 1
82140 Olching
Germany

For devices manufactured and sold in the USA

Printed in the USA

Manufactured for

Hamilton Medical AG
Via Crusch 8, 7402 Bonaduz, Switzerland
☎ +41 58 610 10 20
info@hamilton-medical.com
www.hamilton-medical.com

Distributor in USA

Hamilton Medical, Inc.
4665 Aircenter Circle Reno, NV 89502
☎ (775) 858-3200
Toll-free: (800) 426-6331
info@hamiltonmedical.com

612318/07
2021-11-16

Specifications are subject to change without notice. For all proprietary trademarks and third-party trademarks used by Hamilton Medical AG, see www.hamilton-medical.com/trademarks. © 2021 Hamilton Medical AG. All rights reserved.