
<b>Instructions for use</b>	Nasal cannula, adult/pediatric, single use .....	9
<b>Gebrauchsanweisung</b>	Nasenanüle, Erwachsene/Pädiatrie, für den Einmalgebrauch .....	12
<b>Instructions d'utilisation</b>	Canule nasale, adulte/enfant, à usage unique .....	15
<b>Instrucciones de uso</b>	Cánula nasal, para pacientes adultos y pediátricos, desechable .....	18
<b>Istruzioni per l'uso</b>	Cannula nasale, adulto/pediatrico, monouso .....	21
<b>Instruções de uso</b>	Cânula nasal, adulto/pediátrico, uso único .....	24
<b>使用说明</b>	鼻导管, 成人/儿童, 一次性使用 .....	27
<b>使用説明書</b>	鼻カニューレ、成人/小児用、単回使用 .....	30
<b>Инструкции по эксплуатации</b>	Одноразовая назальная канюля для взрослых/ педиатрических пациентов .....	33
<b>Kullanım Talimatları</b>	Nazal kanül, yetişkin/pediyatrik, tek kullanımlık .....	36

**REF** 10076604, 10076605, 10076606


CE 0197

EAC




**HAMILTON**  
**MEDICAL**  
Intelligent Ventilation since 1983



11



1	Head strap   Kopfband   Bandeau   Banda para la cabeza   Fascia per la testa   Faixa para cabeça   头绑带   ヘッドストラップ   Крепление для головы   Başlık kayışı
2	Adjustment strap   Einstellband   Bandeau ajustable   Banda de ajuste   Cinghia di regolazione   Faixa de ajuste   调节绑带   調節ストラップ   Регулировочный ремень   Ayarlama kayışı   Justeringsstrop   Stelband   Justeringsrem
3	Delivery tube holding clip   Clip zur Befestigung des Schlauchs   Clip de maintien de la tubulure d'administration des gaz   Enganche del tubo de suministro   Clip di fissaggio del tubo di erogazione   Clipe do tubo de fornecimento   输送管路卡夹   送気チューブ保持クリップ   Зажим для крепления подающего шланга   İletim tüpü tutma klipsi
4	Delivery tube   Schlauch   Tubulure d'administration des gaz   Tubo de suministro   Tubo di erogazione   Tubo de fornecimento   输送管路   送気チューブ   Подающий шланг   İletim tüpü
5	Quick release   Schnellentriegelung   Dégagement rapide   Desconexión rápida   Rilascio rapido   Liberação rápida   快速解除扣   クイックリリース   Быстроразъемное соединение   Hızlı serbest bırakma
6	Sealing plug   Verschlussstopfen   Bouchon d'étanchéité   Tapón de sellado   Cappuccio di chiusura   Conector de vedação   密封塞   シーリングプラグ   Уплотнительная заглушка   Sızdırmazlık tıkaçı
7	Face mount, nasal prongs   Gesichtsbefestigungsteil, Nasenprongs   Support facial, embouts nasaux   Montaje facial, adaptadores nasales   Porzione anteriore, nasal prongs   Suporte facial, cateteres nasais   面部支架, 鼻塞   フェイスマウント、鼻ブロング   Лицевое крепление, наконечники назальной каниюли   Yüz bağlantı parçası, nazal pronglar
8	Delivery plug, locking ring   Zufuhranschluss, Sicherungsring   Raccord d'alimentation, bague de verrouillage   Conector de suministro, anillo de bloqueo   Spinotto di erogazione, anello di bloccaggio   Conector de fornecimento, anel de bloqueio   输送塞, 锁环   送気プラグ、ロックリング   Разъем подающего шланга, уплотнительное кольцо   İletim tıkaçı, kilitleme halkası
9	Securing clip   Sicherungsclip   Clip de retenue   Clip de seguridad   Clip di sicurezza   Clipe de segurança   固定夹   固定クリップ   Крепежный зажим   Sabitleme klipsi
10	Connection to inspiratory limb, OD22   Anschluss für den Inspirationsschlauch, AD 22   Connexion à la branche inspiratoire, DE 22   Conexión a la rama inspiratoria, DE22   Collegamento alla branca inspiratoria, D.E. 22   Conexão à alça inspiratória, DE22   连接至吸气肢, 外径 22   吸気回路への接続部、外径22   Соединения для подключения к патрубку вдоха, НД 22   Inspiratuvur parçaya giden bağlantı, DÇ22
11	Adapter, OD15   Adapter, AD 15   Adaptateur, DE 15   Adaptador, DE15   Adattatore, D.E. 15   Adaptador, DE15   接口, 外径 15   アダプタ、外径15   Адаптер, НД 15   Adaptör, DÇ15


4









5

- 
 Partially open nare | Nicht komplett verschlossenes Nasenloch | Narine partiellement ouverte | Fosas nasales parcialmente abiertas | Narice parzialmente aperta | Narina parzialmente aberta | 部分开放鼻孔 | 鼻孔が部分的に開放 | Частично перекрытые ноздри | Kısmen açık burun deliği
- 
 Blocked nare | Blockiertes Nasenloch | Narine bouchée | Fosas nasales bloqueadas | Narice ostruita | Narina bloqueada | 阻塞鼻孔 | 鼻孔が閉塞 | Перекрытые ноздри | Tıkalı burun deliği

10076606 Nasal cannula with adapter, Size S, Adult/Pediatric (orange)

10076605 Nasal cannula with adapter, Size M, Adult/Pediatric (blue)

10076604 Nasal cannula with adapter, Size L, Adult/Pediatric (green)

## Conventions

### ⚠ WARNING

Alerts the user to the possibility of injury, death, or other serious adverse reactions associated with the use or misuse of the device.

### ⚠ CAUTION

*Alerts the user to the possibility of a problem with the device associated with its use or misuse, such as device malfunction, device failure, damage to the device, or damage to other property.*

### NOTICE

Emphasizes information of particular importance.

## 1 Overview

For an overview of the nasal cannula and its components, see the figure on page 2.

### 1.1 Intended use

The nasal cannula is a patient interface intended for the delivery of heated and humidified respiratory gases.

#### NOTICE

- The nasal cannula is intended for use by qualified, trained personnel.
- The intended areas of use are health-care facilities (hospitals and long term acute care hospitals).
- Do *not* use the nasal cannula on patients with occluded or defective nares.
- The nasal cannula is intended for use with patients older than two years.

## 2 Safety information

### ⚠ WARNING

- It is recommended that additional independent monitoring devices, including pulse oximeters measuring SpO<sub>2</sub>, be used at all times. Failure to

monitor the patient may result in loss of therapy, injury, or death.

- Reprocessing of Hamilton Medical single-use products can affect the product properties and may cause injury to the patient. For example, a change to the surface structure during reprocessing may lead to a change in the tear strength or cause actual cracking. Furthermore, an altered surface structure may result in, for example, a microbial aggregation of spores, allergens, and pyrogens, or cause an increase in the number of particles released as a result of chemical changes in the material properties.
- Hamilton Medical does not assume any warranty for the proper function of single-use items if they are reprocessed and reused by the user. Any liability on the part of Hamilton Medical for any damages resulting from the use of reprocessed or reused single-use parts, or use of any non-Hamilton Medical accessories and consumables, is excluded.

- Do *not* use near open flames. *Danger of fire!*
- To prevent loss of therapy, be sure to route the breathing circuit and delivery tube without tension and without any kinks from the device to the patient.
- To prevent disconnection of the delivery tube, only connect to products that comply with ISO 5367.

#### **⚠ CAUTION**

- *Prior to placing the nasal cannula on the patient, turn on the device to ensure that the flow and temperature are appropriate for the patient.*
- *Regularly check for condensation, and drain, if required.*
- *Check the product for damage prior to use. Discard if there is any sign of damage.*
- *Periodically check the position and stability of the prongs.*
- *(USA only) Federal law restricts this device to sale by or on the order of a physician.*

#### **NOTICE**

- Prior to using the nasal cannula, read the high flow therapy device *Instructions for use*.
- Prior to use, check the stability of all connections.

### **3 Setting up and connecting the nasal cannula**

Connecting the nasal cannula comprises the following steps.

Step	See ...
1. Using the sizing guide on the cannula package, select the appropriate size cannula for the patient.	Page 4
(Optional) Connect the delivery tube to the other side of the face mount.	Page 8
2. Connect the nasal cannula to the breathing circuit. Use adapter, if required.	Page 4
3. Turn on the high flow device and check for flow at the prongs.	Page 4
4A. Place the head strap over the patient's head and position the prongs as appropriate.	Page 5
4B. Using the quick release, position the cannula around the patient's head, reconnect the quick release, and position the prongs as appropriate.	Page 5
5. Adjust the head strap. <b>Do not block the nares.</b>	Page 6

#### Step

#### See ...

6. If required, attach the delivery tube to the holder on the head strap. Page 7

7. Connect the securing clip to the patient shirt. Page 7

### **4 Cleaning and replacement**

The cannula is a single-use product that can be used for a maximum of 14 days.

Replace the product in accordance with hospital infection control procedures, or depending on the patient's secretions and nebulization of medication.

Do *not* soak, wash, or sterilize the nasal cannula.

### **5 Delivery**

The components are delivered clean and ready for clinical use.

### **6 Disposal**

Used components must be handled as contaminated. Follow all local, state, and federal regulations with respect to waste management and environmental protection when disposing of used parts.

## 7 Environmental and storage conditions

Environment	Specification
Operating temperature	18°C to 35°C (64°F to 95°F)
Storage temperature	18°C to 35°C (64°F to 95°F)
Transport temperature (for a maximum of 4 weeks)	-10°C to 50°C (14°F to 122°F)
Operating and storage humidity	Less than 95% relative humidity, noncondensing

## 8 Technical specifications and standards

Parameter	Specification
Oxygen concentration	21% to 100%
Maximum flow (l/min)	Size S: 60 Size M: 80 Size L: 100
Flow resistance at 60 l/min (cmH <sub>2</sub> O)	Size S: 17 Size M: 8.5 Size L: 5
Adapter	Connection to OD15*

Parameter	Specification
Connection	OD22*
Maximum operating duration	14 days

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Specifications are subject to change without notice. For all proprietary trademarks (®) and third-party trademarks (§) used by Hamilton Medical AG, see [www.hamilton-medical.com/trademarks](http://www.hamilton-medical.com/trademarks). © 2020 Hamilton Medical AG. All rights reserved.
Hamilton Medical AG
medin Innovations GmbH