Instructions for use  
 Gebrauchsanweisung  
 Instructions d'utilisation  
 Instrucciones de uso

Istruzioni per l'uso  
 Instruções de uso  
 使用说明

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## NIV Full Face Mask vented, non-vented

624942/02 | 2018-04-20








Hamilton Medical AG  
 Via Crusch 8, 7402 Bonaduz, Switzerland  
 +41 58 610 10 20  
 info@hamilton-medical.com  
[www.hamilton-medical.com](http://www.hamilton-medical.com)

**HAMILTON  
 MEDICAL**  
 Intelligent Ventilation since 1983

CE 0197


**HAMILTON  
 MEDICAL**  
 Intelligent Ventilation since 1983

## 9 Pressure/flow characteristic curve

The pressure/flow characteristic curve shows the leak flow as a function of therapy pressure.

Outlet pressure at the therapy device (hPa)

## 10 Materials

Mask component	Material
Mask cushion, forehead cushion	Silicone
Cover for pressure measurement port, spring of forehead support	Silicone
Retaining ring (hard component)	PP (polypropylene)
Retaining ring (soft component)	TPE (thermoplastic elastomer)
Headgear clip	PA (polyamide), POM (polyoxymethylene)
Forehead support, rotating sleeve, mask body, elbow	PA (polyamide), PC (polycarbonate)
Headgear	Elastane, polyester, PU (polyurethane), UBL loop, cotton, PA (polyamide)
Anti-asphyxia valve	TPE, PP (polypropylene)
Quick release cord	POM, polyester, PA (polyamide)

All mask components are free from latex, PVC (polyvinyl chloride), and DEHP (diethylhexyl phthalate).

## 11 Warranty

The manufacturer provides the customer a limited manufacturer warranty on new original products, and on any replacement part fitted by Hamilton Medical AG, in accordance with the warranty conditions applicable to the product in question and in accordance with the warranty periods from date of purchase as listed below. The warranty conditions can be downloaded from [www.hamilton-medical.com](http://www.hamilton-medical.com) on the Internet. We will also send you the warranty conditions on request.

Bear in mind that any claim to warranty and liability shall be void if neither the accessories recommended in the *Instructions for use* nor genuine replacement parts are used. In the event of a claim under warranty, contact Hamilton Medical AG or a distribution partner.

Product	Warranty period
Masks including accessories	6 months

## 12 Declaration of conformity

The product complies with the relevant regulations of Directive 93/42/EEC governing medical devices. The complete text of the declaration of conformity can be obtained from Hamilton Medical AG.

## 13 Ordering information

Description	PN box of 1
<b>Full-face mask, NIV, vented</b>	
Size S	282539
Size M	282538
Size L	282537
<b>Full-face mask, NIV, non-vented</b>	
Size S	282542
Size M	282541
Size L	282540
<b>Other</b>	
Headgear	282632
Quick release cord	282633

## 14 Symbols

Symbol	Description
	Applicable to pediatric/adult patient groups <sup>a</sup>
	Reusable

a. Pediatric patients are understood as children > 30 kg

## 1 Introduction

### 1.1 Intended use

The masks are used for treating sleep apnea and for non-invasive and non-life-sustaining ventilation of patients with respiratory insufficiency (> 30 kg). They serve as a connecting element between the patient and the therapy device. Only use non-vented full face masks in combination with therapy devices equipped with an active exhalation valve.

### 1.2 Contraindications

The mask may not be used, or may be used only with particular caution, if any of the following symptoms are present:

Erosions and ulcerations, skin allergies, flush of the facial skin, pressure marks on the face, claustrophobia, anxiety, deformities of the face or nasopharynx, intake of medication which may induce vomiting, the necessity for immediate intubation, pronounced reflux, or dysfunctional cough reflex. Also observe the contraindications in the *Instructions for use* of your therapy device.

### 1.3 Side effects

Side effects may include: Nasal congestion, dry nose, dry mouth in the morning, feeling of pressure in the sinuses, irritated conjunctiva, skin rashes, pressure marks on the face, irritating noises when breathing. Masks may cause teeth, gums, or jaw to ache.

## 2 Safety

Read these instructions carefully. They are a constituent part of the mask and must be available at all times.

### ⚠ WARNING

- Risk of injury from re-inhaling CO<sub>2</sub>!
- If the mask is used incorrectly, CO<sub>2</sub> may be re-inhaled.
- Vented variant only: Do not close off the exhalation system of the mask.
- Only put on the mask for a prolonged period if the therapy device is running.
- Only use the mask within the quoted therapy pressure range.
- Do not use the mask on patients who are unable to take it off themselves.
- Only use non-vented full face masks in combination with therapy devices equipped with an active exhalation valve.

### ⚠ CAUTION

- Risk of injury from the supply of oxygen!
- Oxygen can become deposited in clothing, bed linen, and hair. In conjunction with smoking, naked flames, or electrical equipment, it can cause fires and explosions.
- Do not smoke.
- Keep away from naked flames.
- Risk of injury if the mask slips! If the mask slips or falls off, the therapy is ineffective.
- Monitor patients with restricted spontaneous respiration.
- Activate low pressure/leakage alarms on the therapy device.

### NOTICE

Risk of injury from anesthetic gases!  
Anesthetic gas may escape through the exhalation system and put third parties at risk. Never use the mask during anesthesia.

## 3 Product description

### 3.1 Overview

Illustrations of the individual components are on the front cover.

#### Key

- 1 Headgear
- 2 Forehead cushion
- 3 Forehead support
- 4 Retaining ring
- 5 Elbow
- 6 Anti-asphyxia valve (vented only)
- 7 Plug for pressure measurement port (vented only)
- 8 Pressure measurement port (vented only)
- 9 Rotating sleeve (vented only)
- 10 Mask body
- 11 Headgear clip
- 12 Mask cushion

### 3.2 Compatible devices

With many combinations of devices, actual pressure in the mask will not correspond to the prescribed therapy pressure, even if the therapy device is displaying the correct pressure. Have the combination of devices adjusted by a physician or specialist dealer so that the actual pressure in the mask corresponds to the therapy pressure.

### 3.3 Exhalation system

#### Vented variants

Vented variants (transparent retaining ring and elbow) have an integrated exhalation system. The retaining ring and mask body are shaped so that there is a gap between these parts. The exhaled air can escape through this gap.

#### Non-vented variants

The non-vented versions (blue retaining ring and elbow) do not have an exhalation system. Only use non-vented full-face masks in combination with therapy devices that have an active exhalation valve and have alarms and safety systems for any device failure.

If you are using the non-vented variant with an external exhalation valve, refer to the associated *Instructions for use*.

### 3.4 Anti-asphyxia valve

#### ⚠ WARNING

- Risk of asphyxia if the anti-asphyxia valve is not working properly!
- Residues may cause the valve to stick and lead to CO<sub>2</sub> being re-inhaled.
- Check before every use that the openings of the anti-asphyxia valve are clear.

If the therapy device fails, the anti-asphyxia valve opens so that the patient can breathe ambient air.

### 3.5 Pressure measurement port

The mask has a pressure measurement port that can be used for measuring therapy pressure or supplying oxygen. If you are not using this port, cover the port with the provided plug.

## 4 Operation

The figures on the inside front and back covers of this guide show you how to adjust, apply, remove, dismantle, and assemble the mask:

Figure 1: Putting on the mask

Figure 2: Adjusting the mask

Figure 3: Removing the mask

Figure 4: Dismantling the mask

Figure 5: Assembling the mask

## 5 Hygienic preparation

#### ⚠ WARNING

- Risk of injury from inadequate cleaning!
- Residues may congest the mask, impair the integrated exhalation system, and jeopardize therapy success.
- For patients with a compromised immune system or particular background of illness, disinfect mask components daily following consultation with the physician.

### 5.1 Cleaning the mask

#### NOTICE

Discolorations of mask components do not impair the functionality of the mask.

- 1 Dismantle the mask (see Figure 4).
- 2 Clean the mask in accordance with the information provided in the following table.

Action	Daily	Weekly
Wash mask components with warm water and a mild detergent.	X	
When washing mask parts, clean thoroughly with a cloth or a soft brush.		X
Wash headgear by hand.		X

- 3 Rinse all parts with clear water.
- 4 Allow all parts to air-dry.
- 5 Perform a visual inspection.
- 6 If necessary: Replace damaged parts.
- 7 Assemble the mask (see Figure 5).

### 5.2 Cleaning, disinfection, and reprocessing

You must clean, disinfect, and sterilize the mask before use on a new patient.

For detailed reprocessing instructions, use the information provided in the *NIV Masks Reprocessing Guide* (PN 624943). The guide is available online at [www.hamilton-medical.com](http://www.hamilton-medical.com).

A printed copy can be sent to you upon request; contact your Hamilton Medical technical representative.

## 6 Disposal

Used masks may be contaminated. Dispose of all parts according to your institution's protocol. Follow all local, state, and federal regulations with respect to environmental protection.

## 7 Troubleshooting

Fault	Cause	Remedy
Facial pressure pain	Mask too tight	Loosen headgear slightly
Draft in the eyes	Mask too loose	Tighten headgear slightly
	Mask does not fit	Use a different size mask
Therapy pressure is not reached	Mask not correctly adjusted	Re-adjust mask (see Figure 2)
	Mask cushion is damaged	Replace mask cushion
	Tube system leaks	Check connection and fit of tubes
	Air is escaping from pressure measurement port (vented only)	Close pressure measurement port with provided cover
Anti-asphyxia valve incorrectly fitted (vented only)		Fit anti-asphyxia valve correctly (see Figure 5)
Anti-asphyxia valve defective (vented only)		Replace anti-asphyxia valve

## 8 Technical data

Product class according to 93/42/EEC	Ila
Dimensions: (W x H x D)	120 mm x 150 mm x 100 mm
Weight:	
• Vented	98 g
• Non vented	98 g
Therapy pressure range	4 hPa - 40 hPa
Hose connection. Cone according to EN ISO 5356-1:	
• Vented	Ø 22 mm (male)
• Non-vented	Ø 22 mm (female)
Pressure measurement port	Ø 4 mm
Width of headgear clips	24 mm
Temperature range:	
• Operation	5°C to 40°C
• Storage	-20°C to 70°C
Flow resistance, vented:	
• at 50 l/min	0.05 hPa
• at 100 l/min	0.28 hPa
Flow resistance, non-vented:	
• at 50 l/min	0.08 hPa
• at 100 l/min	0.32 hPa
Anti-asphyxia valve flow resistance (vented only):	
• Inspiration at 50 l/min	0.5 hPa
• Exhalation at 50 l/min	0.5 hPa
Switching pressure for anti-asphyxia valve (vented only):	
• Open	1 hPa
• Closed	2 hPa
Quoted two-figure noise emission value according to ISO 4871 (vented only):	
• Sound pressure level	19 dB(A)
• Sound power level	27 dB(A)
• Uncertainty factor	3 dB(A)
Service life	Up to 12 months <sup>a</sup>
Standards applied	EN ISO 17510-2: 2009

- a. The materials of the mask deteriorate if exposed to, among other agents, aggressive detergents. In individual cases, therefore, it may be necessary to replace mask parts more frequently.

Hamilton Medical reserves the right to make design modifications.

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